

GLOBAL ANALYTICAL SOLUTIONS



**GERSTEL**

AppNote 2/1992

## Applications of Mass Flow Controlled Multi Column Switching in On-Line Capillary GC/MS

Andreas Hoffmann, Ralf Bremer, Friedhelm Rogies  
*Gerstel GmbH & Co.KG, Eberhard-Gerstel-Platz 1,  
D-45473 Mülheim an der Ruhr, Germany*

Jacques A. Rijks<sup>1,2,3</sup>, Jacek Staniewski<sup>2,3</sup>

<sup>1</sup>*Foundation Education and Development of Chromatography,  
P.O. Box 4990, NL-5604 CD Eindhoven, The Netherlands*

<sup>2</sup>*Poznan Technical University, Institute of Chemical Technology and  
Engineering, Pl. M. Skłodowskiej-Curie 2, 60-965 Poznan, Poland*

<sup>3</sup>*Eindhoven University of Technology, Laboratory of Instrumental  
Analysis, P.O. Box 513, NL-5600 MB Eindhoven, The Netherlands*

### **KEYWORDS**

PTV-injection, multi column switching, cryotrapping,  
capillary GC-MS, industrial applications

## **ABSTRACT**

The separation and analysis of low concentrations of organic compounds in complex sample matrices, such as petroleum products, waste and drinking water, food, beverages and pharmaceutical products is a rather complex analytical problem. Current methods are hampered by insufficient resolution obtained by single capillary columns even if they have rather high plate numbers.

In this paper the potential of a combination of programmed temperature sample introduction and mass flow controlled multi column dual oven capillary gas chromatography and on-line mass spectrometry will be discussed and illustrated.

The effect of cold trapping in between the columns for components with a moderate volatility will be demonstrated for different applications dealing with the determination of trace impurities in various main products such as gasoline, aromatics, steroids and aniline.

## **INTRODUCTION**


Due to ever increasing demands for resolution, sample volume and detection limits in capillary GC nowadays, multi column switching is becoming increasingly important. Essentially these problems are related to the complexity of the samples and the required compatibility of sample size, input band width, column properties, detector specifications and the actual concentrations of the components of interest in the sample.

To realize adequate analysis a number of clean-up and/or enrichment steps, or solvent elimination for diluted samples prior to separation, are needed in order to achieve the required separation efficiency and detection limits. The identification by an on-line coupled MS-system is another breakpoint in the analysis of such complex mixtures, particularly if they elute in low concentrations close together with a major component. Multi column switching is, therefore, an attractive approach to solve these problems, particularly if high and low concentrations have to be determined simultaneously.

The aim of this paper is to demonstrate the possibilities of a combination of temperature programmed sample introduction, multi column dual oven capillary GC, coupled on-line with mass spectrometry in order to enhance identification. The effect of cryotrapping at the inlet of the second column will be emphasized.


## EXPERIMENTAL

*Instrumentation.* The system consists of a combination of a manual or automatic injection (HP 7673, Hewlett-Packard, Avondale, USA), a temperature programmable cold injection system (CIS-3, Gerstel GmbH, Mülheim an der Ruhr, Germany), a multi column switching system (MCS A or MCS P, Gerstel GmbH, Mülheim an der Ruhr, Germany), 2 Ovens (HP 5890 series II, Hewlett-Packard, Avondale, USA) with a cryotrap system (CTS-1, Gerstel GmbH, Mülheim an der Ruhr, Germany) in between and a mass selective detector (HP 5971, Hewlett-Packard, Avondale, USA) as illustrated in **Figure 1**.




**Figure 1.** Schematic diagram of the applied system consisting of a manual or automated injector, a temperature programmable cold injection system with septumless sampling head (2), a GC (3) configured with monitor-FID, column switching device (4) and MCSA or MCSP-pneumatic, connected via a heated transferline with included cryotrap (5) to a second GC (6) with a mass selective detector (7).

The pneumatics of the multi column switching systems used in this study are presented in **Figure 2a** and **b**.



**Figure 2a.** Schematic design of a dual column MCS P system. Solvent flush prior to the first column NV1 and MFC2, both open.




**Figure 2b.** Schematic design of the multi column MCS A system. Solvent flush via splitline, MFC3 open.

The switching pneumatics of both systems are mass flow controlled. The carrier gas flow of the MCS P system is mass flow controlled, the carrier gas flow of the MCS A system is pressure controlled.


## RESULTS

Multi column switching is the most powerful approach for analytical and preparative separation of components which can hardly or not at all be separated in single column systems. This is illustrated below for 5 different applications all dealing with the analysis of trace components, which elute very close to or together with a major component or the determination of trace impurities in highly pure products.

**Example 1.** Multi column determination of benzene and 1-me-cyclopentene in gasoline.



**Figure 3a.** Precolumn chromatogram, 1  $\mu$ l, splitless, marked components transferred to the main column.



**Figure 3b.** Main column chromatogram (TIC).

*Analysis conditions.*

**System:**

MCS A, 2 columns, 2 ovens, CTS, MSD

**Columns:**

Pre-column in GC 1    50 m HP-1 (Hewlett-Packard)  
 $d_i = 0.32$  mm     $d_f = 1.05$   $\mu$ m.  
Main column in GC 2    50 m Wg-11 (WGA)  
 $d_i = 0.25$  mm     $d_f = 0.2$   $\mu$ m.

**Pneumatics:**

Carriergas    He     $p_i = 100$  kPa    split x:50  
Control flow     $p_c = 30$  kPa    10 ml/min  
FID     $H_2$ , 30 ml/min    Air, 300 ml/min  
 $N_2$ , 30 ml/min

**Temperatures:**

CIS	40°C;	$\nearrow$ 320°C;	12°C/s.
Oven 1	50°C;	$\nearrow$ 280°C;	5°C/min.
Oven 2	50°C;	$\nearrow$ 120°C;	5°C/min.
CTS	220°C;	$\searrow$ -150°C;	12°C/s;
		$\nearrow$ 220°C;	12°C/s.


**Detectors:**

Monitor detector    FID  
Main detector    MSD    Scan 10-350 amu


**Compounds:**

1. 1-Me-Cyclopentene
2. Benzene

**Example 2.** Multi column analysis of m-,p-xylene in ethylbenzene.



**Figure 4a.** Precolumn chromatogram, 1  $\mu$ l, split ratio x:30, marked compound flushed.



**Figure 4b.** Main column chromatogram, without CTS.

*Analytical conditions.*

**System:**

MCS A, 2 columns, 2 ovens, with and without CTS

**Columns:**

Pre-column in GC-1	25 m Ultra-1 (Hewlett-Packard)
	$d_i = 0.32$ mm $d_f = 1.02$ $\mu$ m.
Main column in GC 2	30 m CP-Wax (Chrompack)
	$d_i = 0.32$ mm $d_f = 1.2$ $\mu$ m.

**Pneumatics:**

Carriergas	He	$p_i = 130$ kPa	split x:40
Control flow		$p_c = 50$ kPa	10 ml/min
FID	H <sub>2</sub> , 30 ml/min	Air, 300 ml/min	
	N <sub>2</sub> , 30 ml/min		

**Temperatures:**


CIS	60°C;	$\nearrow 260^\circ\text{C}$ ;	12°C/s.
Oven 1	70°C.		
Oven 2	50°C;	$\nearrow 180^\circ\text{C}$ ;	3°C/min.
CTS	200°C;	$\searrow -150^\circ\text{C}$ ;	12°C/s;
		$\nearrow 200^\circ\text{C}$ ;	12°C/s.

**Detectors:**

Monitor detector	FID
Main detector	FID


**Compounds:**

1. Benzene
2. Toluene
3. Ethylbenzene
4. m-Xylene
5. p-Xylene
6. Styrene




**Figure 4c.** Main column chromatogram, with CTS.

### Example 3. Multi column analysis of by-products in steroids.



**Figure 5a.** Single column analysis in split (A) and on-column (B) mode.



**Figure 5b.** Multi column analysis (TIC), main compounds flushed.

#### Analysis conditions.

##### System:

MCS P, 2 columns, 2 ovens, MSD

##### Columns:

Pre-column in GC 1      10 m SB-1 (IAS)  
    $d_i = 0.32$  mm       $d_f = 0.1$   $\mu$ m.  
Main column in GC 2      50 m SB-1 (IAS)  
    $d_i = 0.32$  mm       $d_f = 0.1$   $\mu$ m.

##### Pneumatics:

Carriergas	He	1.0 ml/min	10ml/min
Control flow	PC 1	0 ml/min	20 ml/min
	PC 2	0 ml/min	8 ml/min
FID	H <sub>2</sub>	30 ml/min	Air, 300 ml/min
	N <sub>2</sub>	30 ml/min	

##### Temperatures:

CIS	60°C;	↗ 320°C;	12°C/s.
Oven 1	50°C;	↗ 280°C;	15°C/min.
Oven 2	50°C;	↗ 280°C;	10°C/min.
Transferline	280°C.		


##### Detectors:

Monitor detector	FID	
Main detector	MSD	Scan 30-450 amu


##### Compounds:

1. Steroid X
2. Steroid Y

**Example 4.** Multi column analysis of impurities in aniline.



**Figure 6a.** Precolumn chromatogram, 1  $\mu$ l, split ratio x:10, marked compound flushed.



**Figure 6b.** Main column chromatogram (TIC), without CTS.

*Analytical conditions.*

**System:**

MCS A, 2 columns, 2 ovens, with and without CTS

**Columns:**

Pre-column in GC-1	25 m OV-17 (home made)
	$d_i = 0.32$ mm $d_f = 1.0$ $\mu$ m.
Main column in GC 2	50 m HP-1(Hewlett-Packard)
	$d_i = 0.32$ mm $d_f = 1.05$ $\mu$ m.

**Pneumatics:**

Carriergas	He	$p_i = 60$ kPa	split x:15
Control flow		$p_c = 45$ kPa	10 ml/min
FID	H <sub>2</sub> , 30 ml/min	Air, 300 ml/min	
	N <sub>2</sub> , 30 ml/min		

**Temperatures:**


CIS	60°C;	$\nabla$ 260°C;	12°C/s.
Oven 1	50°C;	$\nabla$ 200°C;	10°C/min.
Oven 2	40°C;	$\nabla$ 260°C;	10°C/min.
CTS	220°C;	$\nabla$ -150°C;	12°C/s;
		$\nabla$ 220°C;	12°C/s.

**Detectors:**

Monitor detector	FID
Main detector	MSD      Scan 10-300 amu

**Compounds:**


1. Benzene
2. Cyclohexane
3. Cyclohexylamine
4. Cyclohexanol
5. Phenol
6. Aniline
7. Toluidine
8. Nitrobenzene
9. Dicyclohexylamine




**Figure 6c.** Main column chromatogram (TIC), with CTS.



### Example 5. Multi column analysis of impurities in styrene.



**Figure 7a.** Precolumn chromatogram, 0.5  $\mu$ l, split ratio x:30, marked compound flushed.



**Figure 7b.** Main column chromatogram, with CTS.

#### Analytical conditions.

##### System:

MCS A, 2 columns, 2 ovens, CTS

##### Columns:

Pre-column in GC-1	25 m Ultra-2 (Hewlett-Packard)
	$d_i = 0.32$ mm $d_f = 0.52$ $\mu$ m.
Main column in GC 2	30 m DB-Wax (Restec)
	$d_i = 0.32$ mm $d_f = 0.25$ $\mu$ m.

##### Pneumatics:

Carrier gas	He	$p_i = 120$ kPa	split x:30
Control flow		$p_c = 65$ kPa	10 ml/min
FID	H <sub>2</sub> , 30 ml/min	Air, 300 ml/min	
	N <sub>2</sub> , 30 ml/min		

##### Temperatures:

CIS	60°C;	$\nearrow 260^\circ\text{C};$	12°C/s.
Oven 1	60°C;	$\nearrow 150^\circ\text{C};$	5°C/min.
Oven 2	40°C;	$\nearrow 70^\circ\text{C};$	5°C/min;
		$\nearrow 180^\circ\text{C};$	10°C/min.
CTS	150°C;	$\searrow -150^\circ\text{C};$	12°C/s;
		$\nearrow 150^\circ\text{C};$	12°C/s.

##### Detectors:

Monitor detector	FID
Main detector	FID

## CONCLUSIONS

Mass flow controlled dual oven multi column switching in combination with on-line coupled MS is a user friendly, very powerful tool for the analysis of complex samples with a wide range of volatilities, polarities and concentrations. High quality mass spectra are obtained after enrichment of volatile trace components.

Cryotrapping at the inlet of the analytical column results in a significant improvement of the separation of critical peak pairs and the reproducibility of the analysis.



### GERSTEL GmbH & Co. KG

Eberhard-Gerstel-Platz 1  
45473 Mülheim an der Ruhr  
Germany

- +49 (0) 208 - 7 65 03-0
- +49 (0) 208 - 7 65 03 33
- gerstel@gerstel.com
- www.gerstel.com

## GERSTEL Worldwide

### GERSTEL, Inc.

701 Digital Drive, Suite J  
Linthicum, MD 21090  
USA

- +1 (410) 247 5885
- +1 (410) 247 5887
- sales@gerstelus.com
- www.gerstelus.com

### GERSTEL AG

Wassergrabe 27  
CH-6210 Sursee  
Switzerland

- +41 (41) 9 21 97 23
- gerstelag@ch.gerstel.com
- www.gerstel.ch

### GERSTEL K.K.

1-3-1 Nakane, Meguro-ku  
Tokyo 152-0031  
SMBC Toritsudai Ekimae Bldg 4F  
Japan

- +81 3 5731 5321
- +81 3 5731 5322
- info@gerstel.co.jp
- www.gerstel.co.jp

### GERSTEL LLP

Level 25, North Tower  
One Raffles Quay  
Singapore 048583

- +65 6622 5486
- +65 6622 5999
- SEA@gerstel.com
- www.gerstel.com

### GERSTEL Brasil

Av. Pascoal da Rocha Falcão, 367  
04785-000 São Paulo - SP Brasil

- +55 (11)5665-8931
- +55 (11)5666-9084
- gerstel-brasil@gerstel.com
- www.gerstel.com.br



Awarded for the  
active pursuit of  
environmental sustainability

Information, descriptions and specifications in this Publication are subject to change without notice. GERSTEL, GRAPHPACK and TWISTER are registered trademarks of GERSTEL GmbH & Co. KG.

© Copyright by GERSTEL GmbH & Co. KG